MP Biomedicals

Case study - Bones

Mouse bone samples grinding with FastPrep-24[™] homogenizer and Lysing Matrix tubes.

Sarah Beck-Cormier, CR CNRS.

INSERM U791 - LIOAD - Laboratoire d'Ingénierie Ostéo-Articulaire et Dentaire. STEP group «Skeletal Tissue Engineering and Physiopathology». Nantes, France.

Overview

- Keywords: Mouse bone samples, RNA extraction, FastPrep-24™ homogenizer, gene expression
- · Aim of the study: High quality RNA extraction from mouse bone samples
- Application: qPCR
- Sample Name: Mouse Calvaria, long bones, chondrocytes, Osteoblasts (OB).
- Material: FastPrep-24™ homogenizer, Lysing Matrix S (metal beads), Lysing Matrix D (ceramic beads)
- Buffer: RNA extraction buffer

Protocol and Parameters

- 1. Put bone samples immediately in liquid nitrogen after sampling
- 2. Place bone sample in Lysing Matrix tube containing 6 metal beads (Lysing Matrix S)
- 3. Add 350 μI of RNA extraction buffer per sample
- 4. Load tube in the FastPrep-24[™] homogenizer and process 2 x 15 sec at speed setting of 6m/s with 5 min intermediate incubation on ice
- 5. Centrifuge Lysing Matrix tube at 12 000 x g, 5 min at 4°C
- 6. Transfer supernatant to a new 2 ml microcentrifuge tube and follow RNA extraction according to RNA extraction protocol

=	1	2	3	4	5
				1	
	-	-	-	-	-

Results

1-4: 250-300 ng RNA from mineralized OB (Lysing Matrix D)

5: 200 ng RNA from mineralized chondrocytes (Lysing Matrix S)

Conclusion

- Bone sample grinding is challenging. This study shows that the use of FastPrep-24[™] homogenizer in combination with metal beads is **highly performant** for this application.
- High RNA quality and yield are extracted from mineralized Osteoblasts and Chondrocytes using FastPrep-24[™], Lysing Matrix D & Lysing Matrix S tubes.
- FastPrep® technology is the most adapted for bone & Forensics studies.

Successful sample preparation using the MP Biomedicals FastPrep® product line has been highlighted in thousands of scientific articles. To access articles and other materials, visit www.mpbio.com/FastPrepLibrary.

www.mpbio.com