MP Biomedicals

Case study - Feces

Detection of Kudoa septempunctata 18S Ribosomal DNA in Patient Fecal Samples from Novel Food-Borne Outbreaks Caused by Consumption of Raw Olive Flounder.

Tetsuya Harada, Takao Kawai, Michio Jinnai, Takahiro Ohnishi, Yoshiko Sugita-Konishi and Yuko Kumeda. J. Clin. Microbiol. 2012. Vol 50.

Overview

- Keywords: Food-borne disease, Parasites identification, Human feces, qPCR, K.septempunctata
- Aim of the study: Identification of a standard method for DNA extraction from fecal parasites
- Application: Quantitative PCR
- Sample name: Human fecal sample
- Material: FastDNA[™] Spin Kit for Soil containing Lysing Matrix E (MP Biomedicals), QIAamp DNA stool minikit (Qiagen), UltraClean fecal DNA kit (Mo Bio)
- Buffer: provided with each of the three commercial DNA extraction kits

Protocol and Parameters

To compare the amount of K.septempunctata (parasites) DNA extracted using the three kits.

- 1. 200 mg of each sample and 200 µl of DNA elution buffer were used during the extraction procedure for each kit.
- 2. Extracted DNA was stored at -20°C until use.

Conclusion

- The FastDNA™ Spin Kit for Soil showed to be the best DNA extraction method providing the highest PCR amplification.
- The FastPrep® technology gives higher yields and increases detection limit threshold of PCR. FastDNA™ Spin Kit for Soil is the most efficient method for extracting parasites DNA from fecal samples.

Successful sample preparation using the MP Biomedicals FastPrep® product line has been highlighted in thousands of scientific articles. To access articles and other materials, visit <u>www.mpbio.com/FastPrepLibrary</u>.

www.mpbio.com