MP Biomedicals

Case study - Seeds

Metabolite profiling and quantitative genetics of natural variation for flavonoids in Arabidopsis.

Jean-Marc Routaboul, Christian Dubos, Gilles Beck, Catherine Marquis, Przemyslaw Bidzinski, Olivier Loudet and Loic Lepiniec. Journal of Experimental Botany. 2012. Vol 63

Overview

- Keywords: Arabidopsis, flavonoids, metabolite profiling, natural variation, quantitative trait loci
- . Aim of the study: Identification of genes controlling flavonoid metabolism
- Application: Flavonoid extraction and analysis by LC-MS
- Sample name: Arabidopsis
- Sample type: Seeds
- Material: FastPrep-24™ homogenizer
- Buffers: Acetonitrile/Water (3/1; v/v) or Methanol/Acetone/water/Trifluoroacetic acid (30/42/28/0.05; v/v/v/v)

Protocol and Parameters

- 1. Three representative seed aliquots from the three biological repeats were pooled before flavonoid extraction.
- 2. All seed samples were ground for 90 s at maximum speed with a FastPrep-24™ homogenizer in 1 ml of solvent mixe.
- 3. A 4 μg aliquot of apigenin was added as an internal standard.
- 4. Following centrifugation, the pellet was extracted further with 1 ml of the same solvent mixes overnight at 4°C.

Conclusion

• The use of FastPrep-24™ system succeeded in full homogenization of seeds, allowing flavonoid extraction, quantification and complete analysis of their metabolism by LC-MS method.

Successful sample preparation using the MP Biomedicals FastPrep® product line has been highlighted in thousands of scientific articles. To access articles and other materials, visit www.mpbio.com/FastPrepLibrary.



